

Product Datasheet

Haptoglobin (HP) antibody (orb1311856)

Description	Mouse monoclonal antibody to Haptoglobin (HP)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, WB
Immunogen	Full length human recombinant protein of human HP (NP_005134) produced in HEK293T cell.
Target	HP
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	43.3 kDa
Uniprot ID	P00738
NCBI	NM_005143
Dilution Range	WB 1:500~2000, IF 1:100
Expiration Date	12 months from date of receipt.