


Product Datasheet

Cytokeratin 19 (KRT19) antibody (orb1311871)

Description	Mouse monoclonal antibody to Cytokeratin 19 (KRT19)
Species/Host	Mouse
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, WB
Immunogen	Recombinant fragment expressed in E.coli corresponding to amino acids 240-390 of human CK19
Target	KRT19
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	43.9 kDa
Uniprot ID	P08727
NCBI	NM_002276
Dilution Range	WB 1:1000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...