

Product Datasheet

SARS-CoV-2 N Protein antibody (orb1311883)

Description	Rabbit polyclonal antibody to SARS-CoV-2 N Protein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Full length recombinant SARS-CoV-2 N protein produced in E.coli.
Target	N Protein
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47kDa
Dilution Range	ELISA 1:5000-1:10000, WB 1:2000, IHC 1:500-1:1000, ICC 1:500-1:1000
Expiration Date	12 months from date of receipt.

