


Product Datasheet

UBE2NL antibody (orb1311886)

Description	Rabbit polyclonal antibody to UBE2NL
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human UBE2NL
Target	UBE2NL
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17.2 kDa
Uniprot ID	Q5JXB2
NCBI	NM_001012989
Dilution Range	WB: 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...