


Product Datasheet

Synaptopodin (SYNPO) antibody (orb1311890)

Description	Rabbit polyclonal antibody to Synaptopodin (SYNPO)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human SYNPO
Target	SYNPO
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	73.5 kDa
Uniprot ID	Q8N3V7
NCBI	NM_001109974
Dilution Range	WB: 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...