


Product Datasheet

SAPS2 (PPP6R2) antibody (orb1311892)

Description	Rabbit polyclonal antibody to SAPS2 (PPP6R2)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human PPP6R2
Target	PPP6R2
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	101.1 kDa
Uniprot ID	O75170
NCBI	NM_014678
Dilution Range	WB: 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...