


Product Datasheet

RMI2 antibody (orb1311894)

Description	Rabbit polyclonal antibody to RMI2
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human RMI2
Target	RMI2
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15.7 kDa
Uniprot ID	Q96E14
NCBI	NM_152308
Dilution Range	WB: 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...