


Product Datasheet

INCENP antibody (orb1311909)

Description	Rabbit polyclonal antibody to INCENP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human INCENP
Target	INCENP
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105.2 kDa
Uniprot ID	Q9NQS7
NCBI	NM_001040694
Dilution Range	WB: 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...