


Product Datasheet

NMDAR2A (GRIN2A) antibody (orb1311916)

Description	Rabbit polyclonal antibody to NMDAR2A (GRIN2A)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human GRIN2A
Target	GRIN2A
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	165.28 kDa
Uniprot ID	Q12879
NCBI	NM_000833
Dilution Range	WB: 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...