


Product Datasheet

AGXT antibody (orb1311939)

Description	Rabbit polyclonal antibody to AGXT
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human AGXT
Target	AGXT
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42.8 kDa
Uniprot ID	P21549
NCBI	NM_000030
Dilution Range	WB: 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...