


Product Datasheet

Aldehyde Oxidase (AOX1) antibody (orb1311969)

Description	Rabbit polyclonal antibody to Aldehyde Oxidase (AOX1)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human AOX1
Target	AOX1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	147.7 kDa
Uniprot ID	Q06278
NCBI	NM_001159
Dilution Range	WB: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...