


Product Datasheet

CD62P (SELP) antibody (orb1311980)

Description	Rabbit polyclonal antibody to CD62P (SELP)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human SELP
Target	SELP
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90.78 kDa
Uniprot ID	P16109
NCBI	NM_003005
Dilution Range	WB: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...