

Product Datasheet

VR1 (TRPV1) antibody (orb1311992)

Description	Rabbit polyclonal antibody to VR1 (TRPV1)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human TRPV1
Target	TRPV1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	94.8 kDa
Uniprot ID	Q8NER1
NCBI	NM_080704
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.

