


Product Datasheet

MAX antibody (orb1312013)

Description	Mouse monoclonal antibody to MAX
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human MAX (NP_002373) produced in E.coli.
Target	MAX
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	18.1 kDa
Uniprot ID	P61244
NCBI	NM_002382
Dilution Range	WB1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of overexpressed I...