

Product Datasheet

GFAP antibody (orb1312041)

Description	Mouse monoclonal antibody to GFAP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human recombinant protein fragment corresponding to amino acids 292-432 of human GFAP (NP_002046) produced in E.coli.
Target	GFAP
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	49.9 kDa
Uniprot ID	P14136
NCBI	NM_002055
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.