


Product Datasheet

SARS-CoV-2 N Protein antibody (orb1312084)

Description	Mouse monoclonal antibody to SARS-CoV-2 N Protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Full length recombinant SARS-CoV-2 N protein produced in HEK293.
Target	SARS-CoV-2 N Protein
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	47kDa
Dilution Range	WB 1:500-1000, ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of full length rec...