

Product Datasheet

S Protein antibody (orb1312097)

Description	Mouse monoclonal antibody to S Protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 Spike-RBD protein produced in E.coli.
Target	S Protein
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	34kDa
NCBI	YP_009724390
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.