

Product Datasheet

NGAL (LCN2) antibody (orb1312135)

Description	Mouse monoclonal antibody to NGAL (LCN2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human recombinant protein fragment corresponding to amino acids 21-198 of human NGAL (NP_005555) produced in E.coli.
Target	LCN2
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	22.6 kDa
Uniprot ID	P80188
NCBI	NM_005564
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.