


## Product Datasheet

### SARS-CoV-2 N Protein antibody (orb1312143)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to SARS-CoV-2 N Protein                |
| <b>Species/Host</b>        | Mouse                                                            |
| <b>Reactivity</b>          | Virus                                                            |
| <b>Conjugation</b>         | Unconjugated                                                     |
| <b>Tested Applications</b> | ELISA, WB                                                        |
| <b>Immunogen</b>           | Recombinant full length SARS-CoV-2 N protein produced in E.coli. |
| <b>Target</b>              | N Protein                                                        |
| <b>Note</b>                | For research use only                                            |
| <b>Isotype</b>             | IgG2a                                                            |
| <b>Clonality</b>           | Monoclonal                                                       |
| <b>Dilution Range</b>      | WB 1:500-1000, ELISA 1:5000-10000                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |



Western blot analysis of recombinant full length SARS-CoV-2 N protein