


Product Datasheet

Desmin (DES) antibody (orb1312220)

Description	Mouse monoclonal antibody to Desmin (DES)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human Desmin (NP_001918). The exact sequence is proprietary.
Target	DES
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	53.5 kDa
Uniprot ID	P17661
NCBI	NM_001927
Dilution Range	IHC 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...