

Product Datasheet

PGP9.5 (UCHL1) antibody (orb1312255)

Description	Mouse monoclonal antibody to PGP9.5 (UCHL1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Full length human recombinant protein of human UCHL1 (NP_004172) produced in E.coli.
Target	UCHL1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	24.8 kDa
Uniprot ID	P09936
NCBI	NM_004181
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.