

Product Datasheet

FABP3 antibody (orb1312342)

Description	Mouse monoclonal antibody to FABP3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Full length human recombinant protein of human FABP3 (NP_004093) produced in E.coli.
Target	FABP3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	14.9 kDa
Uniprot ID	P05413
NCBI	NM_004102
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.