


Product Datasheet

Calponin (CNN1) antibody (orb1312408)

Description	Mouse monoclonal antibody to Calponin (CNN1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 67-257 of human Calponin1 (NP_001290) produced in E.coli.
Target	CNN1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	33.2 kDa
Uniprot ID	P51911
NCBI	NM_001299
Dilution Range	IHC 1:800
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...