


Product Datasheet

CD28 antibody (orb1312477)

Description	Mouse monoclonal antibody to CD28
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human CD28 (NP_006130) produced in HEK293T cell.
Target	CD28
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	25.1 kDa
Uniprot ID	P10747
NCBI	NM_006139
Dilution Range	WB 1:500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...