


Product Datasheet

Nephrin (NPHS1) antibody (orb1312517)

Description	Mouse monoclonal antibody to Nephrin (NPHS1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1077-1241 of human NPHS1 (NP_004637) produced in E.coli.
Target	NPHS1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	134.6 kDa
Uniprot ID	O60500
NCBI	NM_004646
Dilution Range	WB 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...