


Product Datasheet

YAP1 antibody (orb1312552)

Description	Mouse monoclonal antibody to YAP1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 347-504 of human YAP1 (NP_001123617) produced in E.coli.
Target	YAP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	54.3 kDa
NCBI	NM_001130145
Dilution Range	WB 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Equivalent amou...