


Product Datasheet

LARP1 antibody (orb1312556)

Description	Mouse monoclonal antibody to LARP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 315-615 of human LARP1 (NP_291029) produced in E.coli.
Target	LARP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	123.3 kDa
Uniprot ID	Q6PKG0
NCBI	NM_033551
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...