


Product Datasheet

CD47 antibody (orb1312582)

Description	Mouse monoclonal antibody to CD47
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Full length human recombinant protein of human CD47 (NP_001768) produced in HEK293T cell.
Target	CD47
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.21 kDa
Uniprot ID	Q08722
NCBI	NM_001777
Dilution Range	FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
living 293T ...