


Product Datasheet

CTLA4 antibody (orb1312600)

Description	Mouse monoclonal antibody to CTLA4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Full length human recombinant protein of human CTLA4 (NP_005205) produced in HEK293T cell.
Target	CTLA4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	24.66 kDa
Uniprot ID	P16410
NCBI	NM_005214
Dilution Range	FLOW 1:400
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of living 293T ...