


Product Datasheet

Sortilin (SORT1) antibody (orb1312697)

Description	Mouse monoclonal antibody to Sortilin (SORT1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 34-300 of human SORT1 (NP_002950) produced in E.coli.
Target	SORT1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q99523
NCBI	NM_002959
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...