


Product Datasheet

p75 NGF Receptor (NGFR) antibody (orb1312735)

Description	Mouse monoclonal antibody to p75 NGF Receptor (NGFR)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Full length human recombinant protein of human Ngfr (NP_002498) produced in HEK293T cell.
Target	NGFR
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	45.18 kDa
Uniprot ID	P08138
NCBI	NM_002507
Dilution Range	FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
living 293T
...