


Product Datasheet

Cytoglobin (CYGB) antibody (orb1312766)

Description	Mouse monoclonal antibody to Cytoglobin (CYGB)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CYGB (NP_599030) produced in HEK293T cell.
Target	CYGB
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21.2 kDa
Uniprot ID	Q8WWM9
NCBI	NM_134268
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...