


Product Datasheet

C20orf79 (SCP2D1) antibody (orb1312780)

Description	Mouse monoclonal antibody to C20orf79 (SCP2D1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human C20orf79 (NP_848578) produced in E.coli.
Target	SCP2D1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	17.5 kDa
Uniprot ID	Q9UJQ7
NCBI	NM_178483
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...