


Product Datasheet

Cathepsin L (CTSL) antibody (orb1312789)

Description	Mouse monoclonal antibody to Cathepsin L (CTSL)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 114-333 of human CTSL (NP_001903) produced in E.coli.
Target	CTSL
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	37.53 kDa
Uniprot ID	P07711
NCBI	NM_001912
Dilution Range	WB 1:500, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...