


## Product Datasheet

### ZNF226 antibody (orb1312829)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to ZNF226                                                      |
| <b>Species/Host</b>        | Mouse                                                                                    |
| <b>Reactivity</b>          | Human                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | WB                                                                                       |
| <b>Immunogen</b>           | Full length human recombinant protein of human ZNF226 (NP_001027545) produced in E.coli. |
| <b>Target</b>              | ZNF226                                                                                   |
| <b>Note</b>                | For research use only                                                                    |
| <b>Isotype</b>             | IgG1                                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                               |
| <b>MW</b>                  | 91.7 kDa                                                                                 |
| <b>Uniprot ID</b>          | <a href="#">Q9NYT6</a>                                                                   |
| <b>NCBI</b>                | <a href="#">NM_001032373</a>                                                             |
| <b>Dilution Range</b>      | WB 1:500                                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |



Western blot analysis of HEK293T cells w...