


Product Datasheet

GTF2H3 antibody (orb1312835)

Description	Mouse monoclonal antibody to GTF2H3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-241 of human GTF2H3 (NP_001507) produced in E.coli.
Target	GTF2H3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	34.2 kDa
Uniprot ID	Q13889
NCBI	NM_001516
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...