

Product Datasheet

C1orf19 (TSEN15) antibody (orb1312839)

Description Mouse monoclonal antibody to C1orf19 (TSEN15)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human TSEN15 (NP_443197) produced in E.coli.

Target TSEN15

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 18.5 kDa

Uniprot ID [Q8WW01](#)

NCBI [NM_052965](#)

Dilution Range WB 1:250~500

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...