


Product Datasheet

CAD antibody (orb1312851)

Description	Mouse monoclonal antibody to CAD
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1918-2225 of human CAD (NP_004332) produced in E.coli.
Target	CAD
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P27708
NCBI	NM_004341
Dilution Range	WB 1:250~500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...