


Product Datasheet

MTHFR antibody (orb1312879)

Description	Mouse monoclonal antibody to MTHFR
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-230 of human MTHFR (NP_005948) produced in E.coli.
Target	MTHFR
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	74.4 kDa
Uniprot ID	P42898
NCBI	NM_005957
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...