


Product Datasheet

CD272 (BTLA) antibody (orb1312894)

Description	Mouse monoclonal antibody to CD272 (BTLA)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human recombinant protein fragment corresponding to amino acids 31-157+Mouse Fc of human BTLA (NP_861445) produced in HEK293T cell.
Target	BTLA
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	32.6 kDa
Uniprot ID	Q7Z6A9
NCBI	NM_181780
Dilution Range	FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
living 293T
...