


## Product Datasheet

### TIM 3 (HAVCR2) antibody (orb1312908)

|                            |                                                                                             |
|----------------------------|---------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to TIM 3 (HAVCR2)                                                 |
| <b>Species/Host</b>        | Mouse                                                                                       |
| <b>Reactivity</b>          | Human                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                |
| <b>Tested Applications</b> | FC, WB                                                                                      |
| <b>Immunogen</b>           | Full length human recombinant protein of human HAVCR2 (NP_116171) produced in HEK293T cell. |
| <b>Target</b>              | HAVCR2                                                                                      |
| <b>Note</b>                | For research use only                                                                       |
| <b>Isotype</b>             | IgG1                                                                                        |
| <b>Clonality</b>           | Monoclonal                                                                                  |
| <b>MW</b>                  | 33.2 kDa                                                                                    |
| <b>Uniprot ID</b>          | <a href="#">Q8TDQ0</a>                                                                      |
| <b>NCBI</b>                | <a href="#">NM_032782</a>                                                                   |
| <b>Dilution Range</b>      | WB 1:500, FLOW 1:100                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                             |



Western blot analysis of HEK293T cells w...