


Product Datasheet

TIM 3 (HAVCR2) antibody (orb1312914)

Description	Mouse monoclonal antibody to TIM 3 (HAVCR2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human HAVCR2 (NP_116171) produced in HEK293T cell.
Target	HAVCR2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	33.2 kDa
Uniprot ID	Q8TDQ0
NCBI	NM_032782
Dilution Range	WB 1:500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...