

Product Datasheet

LDL Receptor (LDLR) antibody (orb1312955)

Description Mouse monoclonal antibody to LDL Receptor (LDLR)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications FC

Immunogen Full length human recombinant protein of human LDLR (NP_000518) produced in HEK293T cell.

Target LDLR

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 95.38 kDa

Uniprot ID [P01130](#)

NCBI [NM_000527](#)

Dilution Range FLOW 1:100

Expiration Date 12 months from date of receipt.



Flow
cytometric
analysis of
living 293T

...