


Product Datasheet

Hey L (HEYL) antibody (orb1312965)

Description	Mouse monoclonal antibody to Hey L (HEYL)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human HEYL (NP_055386) produced in E.coli.
Target	HEYL
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	34.9 kDa
Uniprot ID	Q9NQ87
NCBI	NM_014571
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...