


Product Datasheet

GPC1 antibody (orb1312985)

Description	Mouse monoclonal antibody to GPC1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 24-530 of human GPC1 (NP_002072) produced in HEK293T cell.
Target	GPC1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	61.5 kDa
Uniprot ID	P35052
NCBI	NM_002081
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...