


Product Datasheet

CD27 antibody (orb1313000)

Description	Mouse monoclonal antibody to CD27
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CD27 (NP_001233) produced in HEK293T cell.
Target	CD27
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	29.16 kDa
Uniprot ID	P26842
NCBI	NM_001242
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...