

Product Datasheet

Mutarotase (GALM) antibody (orb1313026)

Description Mouse monoclonal antibody to Mutarotase (GALM)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human GALM (NP_620156) produced in HEK293T cell.

Target GALM

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 37.6 kDa

Uniprot ID [Q96C23](#)

NCBI [NM_138801](#)

Dilution Range WB 1:500~2000, IHC 1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...