


Product Datasheet

RGS10 antibody (orb1313029)

Description	Mouse monoclonal antibody to RGS10
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human RGS10 (NP_001005339) produced in HEK293T cell.
Target	RGS10
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	21 kDa
Uniprot ID	O43665
NCBI	NM_001005339
Dilution Range	WB 1:2000, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...