


Product Datasheet

TTC7A antibody (orb1313030)

Description	Mouse monoclonal antibody to TTC7A
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-450 of human TTC7A (NP_065191) produced in E.coli.
Target	TTC7A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	96 kDa
Uniprot ID	Q9ULT0
NCBI	NM_020458
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...