


Product Datasheet

BMP6 antibody (orb1313084)

Description	Mouse monoclonal antibody to BMP6
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 21-374 of human BMP6 (NP_001709) produced in E.coli.
Target	BMP6
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	57 kDa
Uniprot ID	P22004
NCBI	NM_001718
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...